

8/1/92

I.) Chloroplast DNA prep protocol:A.) Follow these steps exactly with NO variations:

- 1.) Grow 750ml of a cell wall deficient strain of C. reinhardtii at 25°C in high light to  $\sim 1 \times 10^7$  cells/ml (in TAP medium).
- 2.) Harvest the cells by centrifugation at 3000 g for 5' at 4°C in Coakridge tubes in JAI.
- 3.) Pour off supernatant.
- 4.) Resuspend cells (the pellet) in 20ml TAP gently using a paint brush.
- 5.) Transfer into a 30ml ml Sarstedt tube & centrifuge at 3000 g for 5' at 4°C (in JAI or JA20).
- 6.) Pour off supernatant.
- 7.) Gently resuspend the cell pellet with a paintbrush in 8ml of "A" buffer (0.1 M NaCl, 50mM EDTA, 20mM Tris-HCl, pH 8.0, at 4°C.
- 8.) Transfer the suspension to a 50ml Fisher tube.
- 9.) Add 0.5ml pronase (10mg/ml in 0.01M Na citrate, pH 5.0 pre-digested at 37°C for 2hr made the same day!!) + mix by swirling.
- 10.) Add 2ml 20% SDS, mix, place tubes in a 50°C 140 bath.)
- 11.) After 45' incubation, add 0.5ml more pronase.
- 12.) Again, after 45' add 0.5ml more pronase.
- 13.) 2 1/2 hours later (incubating at 50°C), the cells should appear brownish. Cool on ice & add 20ml phenol. Mix gently by inversion & leave on ice 20'. Transfer suspension to 30ml Sarstedt tube.
- 14.) Centrifuge the phenol-cell mixture at 7000 rpm for 15' at 4°C.
- 15.) The upper aqueous phase should appear to be a light pink color.
- 16.) Transfer the upper aqueous phase carefully, trying to avoid the interface, to a new Sarstedt tube.

- 17.) glass pasteur pipet.  
Add 2 volumes 100% EtOH, leave in  $-20^{\circ}\text{C}$  overnight.
- 18.) Precipitate the DNA by centrifugation at 8000 R for 10' at  $4^{\circ}\text{C}$ .
- 19.) Invert the tubes on a kimwipe & let sit ~5' to drain the EtOH, wipe off the outer edge of the tube with a kimwipe.
- 20.) Dissolve the pellet in 13ml 10mM Tris, pH 7.5.  
@ let tubes rock 2 hrs.
- 21.) Place tubes in a  $65^{\circ}\text{C}$  H<sub>2</sub>O bath, 1hr to completely redissolve the DNA.
- 22.) Spin 5' at  $4^{\circ}\text{C}$  at 5Krpm to remove particulate matter that has not resuspended.
- 23.) Transfer upper aqueous to a new tube.
- 24.) Add 1.6 ml 1x TEN.  
1x TEN: 10mM Tris-HCl  
10mM EDTA  
150 mM NaCl, pH 8.0.
- 25.) Add 2.4 ml EtBr (0.4mg/ml in 10mM Tris, pH 7.5).
- 26.) Add 23 ml NaI saturated in 1x TEN.  
to saturate a 1x solution of TEN  
add 2g NaI  
1ml 1x TEN.
- 27.) Mix by inversion
- 28.) Take refractive index of each tube & adjust to 1.435!! (do not stray over 1.4385).  
@ This is very important!!  
@ otherwise the gradient is not set up correctly & the DNA will not bend correctly.
- 29.) Transfer to quick-seal tubes
- 30.) Weigh & balance to 0.05 g.
- 31.) Seal caps
- 32.) Spin in Vti50 for 44 hrs at 46Krpm,  $20^{\circ}\text{C}$ .  
@ Set brake to OFF!!

- 33.) Pull topmost chloroplast DNA band with an 18g needle.
- 34.) Extract  $EtBr$  with isogenol alcohol
- 35.) Dialyze DNA overnight against 4L TE (50/10), change once + dialyze against 4L TE (10/1).
- 36.) Precipitate DNA with 2 volumes 100% EtOH,  $\frac{1}{10}$  volume 3M NaOAc, pH 5.2 (overnight at  $-20^{\circ}C$ ).
- 37.) Spin 8K, 10',  $4^{\circ}C$ .
- 38.) Drain pellet
- 39.) Rinse with 1ml 70% EtOH
- 40.) Spin 8K, 10',  $4^{\circ}C$ , drain pellet
- 41.) Dry in speed vac
- 42.) Resuspend DNA in 40 to 1mg/ml.