

DNA PREPARATION FROM *CHLAMYDOMONAS REINHARDTII*

1. 1000 ml of cells grown to stationary (thicker cultures work better).
2. Pellet 5K 2-5 mins.
3. Drain well and resuspend in 8.3 ml of H₂O.
4. 11 ml of resuspended cells are pipetted into 50 ml Corning orange cap tube.
Add 20 ml of SDS-EB.

SDS-EB = 2% SDS
400 mM NaCl
40 mM EDTA
100 mM tris-Cl, 8.0

Note: This solution precipitates on storage.

Mix gently by hand.

5. 15 min. 50°C.
6. Add 32g solid CsCl.
7. Dissolve by shaking on rocker.
Note: Can take a long time.
8. Add 2 ml of 10 mg/ml Ethidium bromide. Mix gently.
9. Load into 38.5 ml quick seal tube. Use funnel provided by Beckman. The solution is very viscous and will not go through a needle.
10. Seal according to Beckman instructions. If tubes aren't full, use CsCl solution (1 g dissolved in 1 ml TE--NOTE: Final volume is greater than 1 ml).
11. 45K, VTi 50, 24 h, 15°-20°
12. Remove band with 5 ml syringe and 18 g needle.
13. Extract with CsCl solution saturated i-Butanol.
14. Dialyze against 0.4M NaCl in TE 24 h.
Change dialysis buffer.
Dialyze against TE 24 h.
15. Phenol extract, chloroform extract.
16. Ethanol ppt, using NH₄OAc, and resuspend in 1.5 mls TE buffer.
17. Add 1.5 g CsCl and 100λ 10 mg/ml EtBr.

18. Spin 90 Krpm 15°C 5 hours in TL 100.3 rotor in TL 100 ultracentrifuge (J. Erickson's lab) using polycarbonate tubes.
19. Pipet out genomic DNA band. Extract EtBr with 1-butanol saturated with CsCl solution.
20. Add to final aqueous phase:
 - 2 volumes TE buffer
 - 6 volumes 100% EthanolLet precipitate on ice.
21. Centrifuge and wash pellet with 70% EtOH.
22. Drain and dry final pellet and resuspend in SH₂O (sterile H₂O).