

FORMALDEHYDE GEL ELECTROPHORESIS OF RNA & RNA Blot Hybridization

see Maniatis (1982), p.202 for reference and see Southern protocol in this manual

Gel Preparation: 3 g + 215 ml H₂O; microwave; stir

In fume hood: While stirring, add 25 ml 10 X MOPS (rt)
10 ml formaldehyde (freshly deionized with AG 501 X 8(D) resin)
Stir to mix well and pour into mold which has comb slightly raised
from normal position. Transfer to gel box before pulling out the
comb.

Sample Preparation: Add appropriate amount (3 – 10 µg/lane) RNA to solution with final concentrations of 50% ultrapure formamide, 6% formaldehyde (deionized) 1 X MOPS. Heat 5 min at 65°C. Add dye.

10 X MOPS Buffer:	0.2 M MOPS	41.8g free acid)Recipe
For 1000 mL	50 mM NaOAc	4.1g)from gene
	10 mM EDTA	20 ml of 0.5 M stock)screen booklet

-titrate to pH 7.0 with NaOH; store at 4°C.

5 x Dye Solution	25% glycerol	10X = 50%
	0.2% bromphenol blue	0.4%
	0.2% xylene cyanol	0.4%

3 µl sample volume:	add in order shown	
	5.1 µl RNA	or 8.1 µl RNA
	3 µl 10 x MOPS	3 µl 10 X MOPS
	12 µl formamide	12 µl formamide
	3.9 µl formaldehyde	3.9 µl formaldehyde
	Heat 5', 65°C	Heat 5', 65°C
	6 µl 5 x Dye Solution	3 µl 10X dye

(If your RNA is concentrated you can make a 20 µl sample by scaling down).

Run gel 4 – 5 hours at 100 V or ON at 40-50 V in 1 X MOPS

Stain

- Stain gel in 500 ml of H₂O containing 50 µl of Ethidium bromide (10 mg/ml stock) 2 h, with gentle shaking
 - Destain in H₂O as necessary to take a picture (usually ON is convenient)
- NOTE: stained gels do not hybridize very well. Therefore, run duplicates, or do not stain.

Blot

Soak gel in 10 X SSC 2 x 15 min on reciprocal shaker

Blot overnight in 10X SSC; use 2 pieces of Whatman 3 MM paper on bottom of gel as wicks (plus 1-2 more pieces to size of gel) – see Southern procedure.

Wet gene screen or other membrane in H₂O and then soak in 10 x SSC for 5'. It is important to remove gel gap from the membrane.

Irradiation

The filter is marked with a felt pen and washed vigorously in 6 X SSC. Place the filter (still wet) on a glass plate and cover it with seran wrap.

Irradiate the filter with 1.6 kJ/m² (1200 μW/cm²) (with a relatively new trans-illuminator at 260 nm this is approximately 5 minute exposure at a distance of 15 cm., but check date of most recent calibration on the illuminator).

Irradiate through Saran wrap – not through glass.

Rinse the filter in distilled water 2 -3 times.

Hybridization and wash procedure from Church & Gilbert (PNAS **81**; 1991-1995, (1984))

Hybridization

Hybridization buffer: 7% SDS

0.5 M NaPO₄, pH 7.2

1% Bovine Serum Albumin (Sigma, fraction V)

Note: 1 M NaPO₄, pH 7.2 is made with 134 grams of Na₂HPO₄ (7 H₂O) plus 4 ml of 85% H₃PO₄ per liter. This is stock solution for hybridization and washing solutions.

Note: This recipe is what I use, but it does not give IM NaP (KH 4/17/91).

Rehybridize the filter for 10 minutes or longer in approximately 0.1 ml hybridization buffer/cm² at 65°C.

Note: I prehybridize and hybridize at 68° (KH 4/17/91)

Remove buffer and add .05 ml/cm² of fresh buffer with probe. Remember to denature probe!!

Hybridize 12 – 16 hours at 65°C.

Wash blot – 6 x/15'/65°C/shaking/100-200 ml of 1% SDS in 40 mM NaPO₄, pH 7.2.

Can wash several blots together without cross contamination. (After hybridization, rinse each filter briefly in a fresh aliquot of 200 ml wash solution (rt); then pool filters in 200 ml for 6 x washes as above).

See southern protocol for exposure.