

OLIGOLABELLING

Feinberg, A.P., Vogelstein, B. (1984) Anal. Biochem. **137**: 266-27

Feinberg, A.P., Vogelstein, B. (1983) Anal. Biochem. **132**: 6-13

Combine in O-ring sealed microfuge tube: 50 – 100 ngs linearized DNA
1.5 μ l Random primers 0.2U/ μ l (Pharmacia)
H₂O to 5 μ l

Cap tightly; 100°C for 3'; plunge into ice.
Microfuge 1''

Add the following to tube:

- 6 μ l 5 x buffer
- 2.5 μ l 5 mg/ml BSA
- 5 μ l 250 μ M dGAT TPs
- 5 μ l α ³²P-dCTP (400-800 Ci/mmol or up to 3000 Ci/mmol)
- 6 μ l H₂O
- 3 – 4 units Klenow (\$0.30/unit)

Incubate @ RT° > 1 h. (Klenow has essentially no activity after a couple hours of incubation.)

Remove 0.5 μ l RXN to 100 μ l TE/0.5 mg/ml carrier DNA to estimate incorporation. See protocol.

SPOT 10 μ l: ~ 170,000
FILTER 10 μ l: ~ 120,000

Can get 8 x 10⁸ cpm/ μ g

	<u>Stock</u>	<u>1 ml</u>
5 x buffer =		
1 M Hepes, 6.6	2 M	500 μ l
250 mM Tris-Cl, 8	1 M	250 μ l
25 mM MgCl ₂	1 M	25 μ l
50 mM 2-ME	14.4 M	3.5 μ l
H ₂ O		221.5

I often use the probe without removing unincorporated nucleotides since I often get 70-80% incorporation.

If incorporation drops, make new buffer.

Boil the probe 10 minutes before using in O-ring tube.

Quench on ice.

Pipet it into hybridization solution.

NOTE: High specific activity nucleotides (3000 Ci/mmol) and probes are more susceptible to radiochemical decomposition. Therefore, use this material only when you need the extra sensitivity.