PROTOCOLE DE CONGELATION A L'AZOTE LIQUIDE

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1) The day before D day, prepare the Nalgene box

Put 250 ml of isopropanol in the box and put the box in the cold room or cold box at 4 C.

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2) On D day : Prepare the cells that you want to freeze : one Culture A per strain

(The quantity for 9 frozen tubes for each strain is given)

- Prepare the cells in a good physiological state 2 streaks per TAP plate. The streaks can be 4 days old, fat and happy cells.
- On D day, resuspend to a high concentration (nice green, about 10^7 per ml) in 10 ml of TAP (this is preculture A)
- Shake 30 min or more so that cells are not aggregated or clumpy and well-dissociated.
- Count cells and calculate the concentration
- Dilute to obtain 10 ml of culture at 6×10^6 per ml in TAP = this is culture A
- Leave Culture A shaking while you are preparing Solution B

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3) Preparation of the solution of methanol : Solution B

(quantity is for 18 frozen tubes, which corresponds to the number of spaces in the Nalgene box)

- Prepare 15 ml of TAP with 10% methanol
- Sterile filter using 4 micron filter
- Distribute 750 μ l solution B in 18 sterile 2 ml cryotubes

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4) Mix of cultures A with solution B in a sterile hood under dim light

- Label the cryotube with strain name, date, etc.
- Take Nalgene box from cold and make sure that isopropanol has not touched in the interior compartment.

Mix A and B : transfer 750 ul of culture A in a tube from step 3 while holding it in the curved palm of your hand, invert the tube 2 times gently.

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6) Second freezing at liquid nitrogen temperature. This has to be done quickly. That is, don’t allow samples from Step 5 to thaw.

The best way to do this is to take the square box that you are using for liquid nitrogen storage and put it in a styrofoam tub with liquid nitrogen. Remove the tubes from the Nalgene box and put the tubes into the long term storage plastic box in the styrofoam container. Make sure that the size of the tubes is compatible with quick removal from one box and placement into the other.

Immediately put the box into the liquid nitrogen tank.

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7) Thawing and Evans Blue staining

After one or two days (or more !) after the day of freezing, thaw the tubes following the instructions on the next page after checking the survival and sterility.

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