

## SOUTHERN BLOT

(Church, G.M., Gilbert, W. *P.N.A.S.* **81**: 1991-1995, April 1984)

1. Run restricted DNA on agarose gel as usual. 1  $\mu\text{g}$  of genomic DNA ( $10^8$  bp genome) will give good signal on overnight exposure with screen if probe is  $10^8$  cpm/ $\mu\text{g}$
2. Denature/Neutralize:
  - .25 M HCl, 15 minutes (for high molecular weight DNA)
  - .5 N NaOH, 1.5 M NaCl; 2 x 20 minutes
  - .5 M Tris-HCl, pH 7.6, 1.5M NaCl; 2 x 20 minutes  
20X SSC, 10 minutes

Soak gene screen in water for 5 minutes, followed by 5 minutes in 10x SSC.

3. Blot:

Overnight in 10 or 20x SSC

Use 2 sheets of 3 mm paper as a wick

Cut additional 3 mm sheets into rectangles the size of the gel

Place one pre-wet sheet on wick

Place gel on Whatman 3 mm rectangle

Place strips of Parafilm on all 4 sides of gel

Place pre-wet membrane on gel

Remove air bubbles. You can use 1 ml plastic pipette as a "squeegee"

Place one more pre-wet 3 mm sheet on membrane

Put a stack of cut 3 mm (1") sheets or cut brown paper towels on top

Place a glass plate over all this

Cover apparatus with Saran wrap

Put Sigma catalog on top as a weight

See p. 384-385 in Maniatis (old one)

Fig. 11-1 (top)

(except there's no need to invert the gel)

4. Irradiation:

The filter is marked with a felt pen and gently immersed in 6x SSC to wash off any gel crap.

Then shake ~20 minutes to reduce future background. Place the filter (still wet) on a glass plate and cover it with Saran wrap. Irradiate the filter with 1.6kJ/m<sup>2</sup> (with a relatively new trans-illuminator at 260 nm, this is approximately a 2 – 2.5 minute exposure at a distance of 15 cm).

Rinse the filter in distilled water 2 – 3 times.

Irradiation substitutes for baking.

*Remember* Irradiate DNA side

and Irradiate through Saran, **not glass**

