

Culture Conditions for Chlamydomonas reinhardtii

A. Culture Media

Stock Solutions (molarities indicated are the final concentrations in tris-acetate-phosphate medium)

Beijerinck's Solution: *6-13-63, Am 4*

NH ₄ Cl	<i>8g</i>	8 g
CaCl ₂ ·2H ₂ O	<i>1ml</i>	1 g
MgSO ₄ ·7H ₂ O	<i>0.82g</i>	2 g
		water to 1000 ml

*7.5 x 10⁻³ M dissolved CaCl₂ in 300ml
 3.4 x 10⁻⁴ M Mg-2Mg → add to 500ml
 4.0 x 10⁻⁴ M Mn → what + MgSO₄ make up to 1l*

Phosphate Buffer:

Na₂HPO₄ (anhydrous)	11.62 g	6.0 x 10⁻⁴ M
KH ₂ PO ₄	7.26 g	4.5 x 10 ⁻⁴ M

or water to 1000 ml

*0.136 M stock K-Pi, pH 7.1
 dilute to 6.33 ml/litre
 pH w/ 'Co Free' KOH*

Tris-Acetate Solution:

Trizma Base (Sigma)	121 g	2.0 x 10 ⁻² M
Glacial acetic acid	50 ml	1.7 x 10 ⁻² M

water to 500 ml

Hutner's Trace Elements Solutions (J. Bacteriol., 52, 213 (1946))

Ethylenediamine tetraacetic acid	50.0 g
ZnSO ₄ ·7H ₂ O	22.0 g
H ₃ BO ₃	11.4 g
MnCl ₂ ·4H ₂ O	5.1 g
FeSO ₄ ·7H ₂ O	5.0 g
CoCl ₂ ·6H ₂ O	1.6 g
CuSO ₄ ·5H ₂ O	1.6 g
(NH ₄) ₆ Mo ₇ O ₂₄ ·4H ₂ O	1.1 g

*50.0 g = 200 ml the separately & add to
 22.0 g = 100 ml of 200 ml of the salts. Do not add
 salt to EDTA. Add EDTA to salts.
 Do not let temp fall below 75°C*

Boil in 750 ml distilled water, cool slightly and bring to pH 6.5-6.8 with KOH (do not use NaOH). The clear solution is diluted to 1000 ml with distilled water and should have a green color which changes to purple on standing. Stable for at least one year.

4-3-63 stir on plate from time to time or continuously
 Media (using stock solutions described above)

Tris-Acetate-Phosphate Medium (pH 7.3):

Beijerinck's	<i>18</i>	<i>32</i>	150 ml
Phosphate	<i>8.33</i>		25 ml
Tracc	<i>1</i>		3 ml
Tris-acetate	<i>10</i>		30 ml
H ₂ O	<i>931</i>	<i>2.792</i>	
			water to 3000 ml

*42. 500
 200 42. 4
 82.3 669 ml.
 4 50 ml
 40. 500 ml
 3.722*

*600ml
 50 ✓
 5
 6
 6
 558 ml*

CHLAMYDOMONAS CULTURE MEDIA

A. SOLUTIONS

1. Phosphate Buffer (2x)

K_2HPO_4 14.34 g
 KH_2PO_4 7.26 g
distilled water 1000 ml

2. Beijerinck's Solution (2x)

NH_4Cl 8 g
 $CaCl_2 \cdot 2H_2O$ 1 g
 $MgSO_4 \cdot 7H_2O$ 2 g

3. Trace Elements (after Hutner, S. H., 1950, Proc. Amer. Phil. Soc. 94: 152-170)

Ethylene diamine (dinitrilo) tetra-acetic acid disodium salt

50.0 g

$ZnSO_4 \cdot 7H_2O$ 22.0 g 16.73 gk. $ZnAc_2$ FW 219

H_3BO_3 11.4 g

$MnCl_2 \cdot 4H_2O$ 5.06g

$FeSO_4 \cdot 7H_2O$ 4.99g 3.54 gcl. $FeCl_2$ FW 199.

$CoCl_2 \cdot 6H_2O$ 1.61g

$CuSO_4 \cdot 5H_2O$ 1.57g 1.26 gAc. $CuAc_2$ FW 200.

$(NH_4)_6Mo_7O_{24} \cdot 4H_2O$ 1.10g

Directions: Add all but the EDTA to 550 ml distilled water, dissolving them one at a time. Put EDTA into 250 ml distilled water, heat until completely dissolved. Bring the first beaker to 70°C, add the EDTA solution. Keeping the total solution at or above 70°C, bring to pH 6.5-6.8 by adding 20% KOH solution (put 20 gm KOH in 100 ml H₂O). (Do not use NaOH, or KOH pellets undissolved. Do not make more alkaline -- it cannot be brought back with acid). Dilute to 1000 ml. Let stand in flask with cotton top until the solution turns purple (at room temperature). Filter out red-brown precipitate. Keep in refrigerator.

B. MEDIA

1. TAP

Trizma 121	2.42 g
glacial acetic acid	1.0 ml
1 M KH ₂ PO ₄ , pH 7.0	1.0 ml
Beijerinck's solution (2x)	50 ml
Trace solution	1 ml
distilled water	950 ml

2. MINIMAL

Beijerinck's solution (2x)	50 ml	plates mm	M4A
Phosphate Buffer (2x)	50 ml	35	35
Trace Solution	1 ml	35	35
Distilled water	900 ml	.7	.7
		630	630
		<u>700ml.</u>	2.3g NaAc
		10.5g agar.	2.8g YE

3. SUPPLEMENTED MEDIA

- a. acetate: minimal plus 2g NaAcetate.3H₂O per litre
- b. yeast acetate: minimal plus ^{3.3}2g Sodium acetate and 4g yeast extract

4. NOTES

- a. for liquid media use as given above. We routinely use 300 ml in a 500 ml flask.
- b. for plates use 1.5% agar, or 15g per litre. Allow 30 ml per plate.
- c. for agar slants use 2% agar, 20g per litre.

C. CARE OF CULTURES

- 1. All stocks are kept in yeast acetate agar slants with screw tops. After inoculation, let them grow up for a week at 500 ft-candles (less for photosynthetic mutants) with the tops loose, then tighten tops and keep in dim light.

2. Mutants currently in use are also kept on several Petrie dishes, yeast extract for wild type and most mutants, TAP for photosynthetic mutants. The plates are transferred every week, or every 3-5 days for photosynthetic mutants. Light is 500 ft-candles for most, 50 or so for the photosynthetics.
3. For biochemical work, flasks of liquid media (usually TAP) are inoculated from plates (the cells first being suspended in a small quantity of sterile water). The shaker has a light level of about 500 ft-candles; for photosynthetic mutants a paper towel around the flask reduces the light by at least 50%.
4. Photosynthetic mutants have a pronounced tendency to revert or become suppressed. They should be cloned and tested on minimal medium frequently. If necessary they can usually be recovered after a backcross to wild type.

6-13-88 Solns

A. Beijerinck's

8g NH_4Cl
1ml CaCl_2 (1g/ml stock)
1g $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$
to 1L

B. Phosphate Buf.

9.25g KH_2PO_4 Aldrich 99.999%
to 500 ml H_2O 22,980-6
= 0.136M stock

8-1 SO_4 , -Cu Beijerinck's

Gold Label [4g NH_4Cl
0.5ml CaCl_2
0.82 $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ (= $8.2 \times 10^{-3}\text{M}$)
to 500 ml w/ H_2O