

August 31<sup>st</sup>, 2006 (Scott Hsieh)

## **Trypsin In-Gel Digest Protocol**

### **Preparation of reagents**

#### 100 mM NH<sub>4</sub>HCO<sub>3</sub>, pH 9.2

Simply dissolve 0.079 g of NH<sub>4</sub>HCO<sub>3</sub> in DI water to a volume of 10 mL. The pH should be around 9.2 after doing this so adjustment is not needed.

#### 10 mM DTT in NH<sub>4</sub>HCO<sub>3</sub> (1.543 mg/mL)

Because of the small amounts of DTT needed to make 10 mM DTT solution it is recommended to make a 1 M DTT stock, which can then be diluted to the needed concentration. A suggested way to do this is to tare an empty microfuge tube, add a small amount of DTT, calculate an appropriate volume of NH<sub>4</sub>HCO<sub>3</sub> to add for the weight of that small amount, then resuspend the DTT in that volume. Make this solution fresh each time you use it.

#### 50 mM Iodoacetamide (IAA) in NH<sub>4</sub>HCO<sub>3</sub> (9.245 mg/mL)

It is recommended to make a 0.5 M IAA (10X or 92 mg/ml) stock, which is then diluted to the needed concentration. You should make this solution fresh each time you need it.

#### CH<sub>3</sub>CN

I usually use CH<sub>3</sub>CN from EMD (AX0151-1)

#### 2% Trifluoroacetic Acid (TFA)

I normally use TFA from Pierce (28903)

10 ng/μL Promega sequencing grade, modified trypsin (Cat no. V5113)

**NOTE:** Volumes given in this protocol are APPROXIMATE. The size of the gel slice will determine how much of each reagent is actually needed. A good rule of thumb is to simply make sure that the solution volume covers at least 75% of the area of the gel. Also it is recommended to work in **500 μL** tubes. Also for steps that require using a speedvac, **do NOT use the speedvac in 5043.** (Acid in samples will ruin pump.) You may use the speedvac in the common area near the - 80°C freezer.

### **Washing**

- 1) If applicable, transfer gel slices into 500 μL tubes before starting. This should reduce the volumes of reagents needed for this protocol.
- 2) If gel slices are dried, allow them to swell for approximately 10 mins in 30-50 μL of 100 mM NH<sub>4</sub>HCO<sub>3</sub>.
- 3) Soak gel slices in 100 μL of a 50:50 mix of 100 mM NH<sub>4</sub>HCO<sub>3</sub>/CH<sub>3</sub>CN. Let stand for 10 mins.

4) Discard supernatant and add 30  $\mu\text{L}$   $\text{CH}_3\text{CN}$ . Let stand for 10 mins. Discard supernatant.

5) Repeat steps 2 and 3 twice or until gel is destained.

6) Dry the gel slice in speedvac for 5-10 mins. **Do NOT use the speedvac in 5043.**

#### **Reduction of disulfides and blocking of free thiols**

1) Add 20  $\mu\text{L}$  of 10 mM DTT to gel slice. Incubate at 60°C for 1 hour on a heating block.

2) Allow samples to cool to room temperature and discard supernatant.

3) Add 20  $\mu\text{L}$  of 50 mM iodoacetamide. Allow samples to incubate in the dark for 45 mins at 45 °C. Discard supernatant.

**Note:** If using volumes different from those given, DTT and iodoacetamide should be added in equal volumes.

#### **Wash gel slices**

1) Add 50  $\mu\text{L}$  of 100 mM  $\text{NH}_4\text{HCO}_3$  and let stand for 10 mins.

2) Discard supernatant and add 50  $\mu\text{L}$   $\text{CH}_3\text{CN}$ . Let stand for 10 mins then discard supernatant.

3) Repeat steps 1 and 2.

4) Dry the samples in the speedvac for 5-10 mins. **Do NOT use the speedvac in 5043.**

#### **Trypsin digest**

1) Swell the gel slices in 10  $\mu\text{L}$  trypsin (refer to step above) and let sit in an ice bath for 45 mins. (Do not worry about covering the gel slice in this step)

2) Add 10-20  $\mu\text{L}$  of 100 mM  $\text{NH}_4\text{HCO}_3$ , and digest overnight at 37 °C. (Do not allow digestion time to exceed 24 hrs.)

3) Add 20  $\mu\text{L}$  of  $\text{H}_2\text{O}$  to each sample and let stand for 5 mins. Transfer supernatant to a new tube.

4) Add 30  $\mu\text{L}$  of 50%  $\text{CH}_3\text{CN}$ /1% TFA to each gel and let stand for 10 mins. Add this new supernatant to the supernatant previously collected.

5) Repeat step 4 two more times.

6) Dry the extracts in the speedvac. (This may take several hours depending on your sample, but be sure not to over dry your extracts.) **Do NOT use the speedvac in 5043.**

7) Store protein at 4°C either dry or resuspend in appropriate buffer dependent on what type of analysis you want to do.

**Note:** It is recommended that you initially resuspend pellets from step 7 in a buffer containing 50% CH<sub>3</sub>CN/ 1% TFA and dry down again to further clean the sample. Then resuspend in buffer that is appropriate for analysis.

For LC-MS/MS: 0.1% Formic acid

For Nanospray MS/MS: 50% CH<sub>3</sub>CN/ 2% Formic acid

For MALDI: 50% CH<sub>3</sub>CN/ 1% TFA